

## Kontrollfacit PNV

ANALYS	Kontroll	-3SD	-2SD	MV	+2SD	+3SD	Lot
<b>Afinion/MicrosEmi</b>							
<b>CRPN</b> mg/L	PNV K4	9	12	16	21	23	K4 1909884, 2103819
	PNV K6	50	57	71	85	92	K6 1909885, 2103820
<b>HemoCue</b>							
<b>Hb</b> g/L	PNV HGL	74	76	79	82	84	HGL 14401, 21801
	PNV HGN	112	114	119	124	127	HGN 14402, 23102
<b>Accu-Chek inform II</b>							
<b>Glukos plasma</b> mmol/L	PNV CG1	1,8	2,0	2,4	2,8	3,0	CG1 20103504, 20100332
	PNV CG2	14,9	15,4	16,5	17,6	18,1	CG2 20103504, 20100332
<b>Clinitek Status</b>							
<b>Glukos</b>	PNV UA1			0			UA1 2106072, 2212108 UA2 2106073, 2301170
	PNV UA2		1		3		
<b>Ketoner</b>	PNV UA1			0			
	PNV UA2		1		4		
<b>Leukocyter</b>	PNV UA1			0			
	PNV UA2		1		4		
<b>Nitrit</b>	PNV UA1			0			
	PNV UA2			1			
<b>Protein</b>	PNV UA1			0			
	PNV UA2		1		3		
<b>Erythrocyter</b>	PNV UA1			0			
	PNV UA2		1		3		
<b>Förklaring:</b>	K4 Seronorm Liquid CRP-kontroll L-1 från Sero K6 Seronorm Liquid CRP-kontroll L-3 från Sero HGL, HGN Hemotrol från Eurotrol CG1, CG2 Accu-Chek Performa Control från Roche UA1, UA2 Liquid QC DipStrip Urinalysis Control från Triolab						